

## Stand Alone Storage Reservoirs

Easy-to-use reservoir specially designed for the storage of purified water. Fully wall-mountable to save valuable bench space.

### Standard Features and Benefits

- Inert polyethylene construction minimizes the risk of organic leachables
- Fully opaque to avoid algal growth
- Integrated 5-level switch control to allow automatic fill-up of reservoir by lab water system
- Sloping base allows complete drainage of water from reservoir
- Composite vent filter eliminates contamination of water from CO<sub>2</sub>, volatile organic compounds or airborne bacteria
- Available in 3 convenient sizes — 25 liter, 40 liter and 75 liter
- Designed for use with PURELAB® systems



### RESERVOIR SPECIFICATIONS

MODEL	HxWxD (in)	Weight (lbs)	Useable Volume
25 Liter Reservoir	16x14x11	11 lbs Empty/ 66 lbs Full	25 Liters
40 Liter Reservoir	21x16x12	18 lbs Empty/ 106 lbs Full	40 Liters
75 Liter Reservoir	27x20x14	26 lbs Empty/ 192 lbs Full	75 Liters

### ORDERING INFORMATION

MODEL	Catalog Number
	<b>For use with MEDICA, Prima and Option</b>
25 Liter Storage Vessel	LA611 (for 7 or 15 lph systems)
40 Liter Storage Vessel	LA612 (for 7 or 15 lph systems)
75 Liter Storage Vessel	LA613 (for 7 or 15 lph systems)
Reservoir 40 Liter	LA587 (for 30,60 or 120 lph systems)
Reservoir 75 Liter	LA590 (for 30,60 or 120 lph systems)
Composite Vent Filter (supplied)	LC136M2
Wall Mounting Bracket (25/40 Liters)	LA591USA
Wall Mounting Bracket (75 Liters)	LA592USA

